



SAFE ANIMAL HANDLING OPERATIONS

ACCORDING TO THE US NATIONAL INSTITUTE FOR OCCUPATIONAL SAFETY AND HEALTH (NIOSH):

- Allergic reactions are among the most common conditions that adversely affect the health of personnel involved in the care and use of animals in research.
- Symptoms of allergies include sneezing, nasal congestion, itchy eyes and cough.
- Some personnel may even develop asthma or other respiratory disorders.
- Personnel should protect themselves from exposure to animals and animal products.
- Seek medical attention in case allergy symptoms develop.

ANIMAL HANDLERS SHOULD TAKE STEPS TO PROTECT THEMSELVES FROM EXPOSURE TO ANIMALS AND ANIMAL PRODUCTS:

- Perform animal manipulations within ventilated workstations or safety cabinets when possible.
- Avoid wearing street clothes while working with animals.
- Leave work clothes at the workplace to avoid potential exposure problems for family members.
- Keep cages and animal areas clean.
- Reduce skin contact with animal products such as dander, serum, and urine by using gloves, lab coats, and approved particulate respirators with faceshields.

CAGE-CHANGING RECOMMENDATIONS FOR SMALL ANIMALS

- Utilize an Animal Containment Workstation which protects both animals and personnel.
- Plan the work session and prepare all items required before starting work.
- Turn on the fan and wipe down the work surface with suitable disinfectant.
- Remove cages from the ventilated rack(s) and change them one at a time in the workstation.
- Spray your gloved hands with disinfectant before removing the rodent from a dirty cage.
- Remove the rodent via holding the mid-base of its tail, and transfer it to a clean cage. If there is more than one rodent in the cage, transfer the rodents one at a time.
- Handle animals with care, as inappropriate handling may provoke animal response and result in possible injury and discomfort.
- Ensure the feeder and water bottle are properly secured.
- Return the clean cage to the correct slot in the ventilated rack system.
- After all cages are changed, perform a final animal count and verify all rodent cages are properly secured in the ventilated rack(s).
- Wipe the work surface with suitable disinfectant and turn off the fan.
- Log all activities as required by your facility.



BEDDING DISPOSAL RECOMMENDATIONS

- If possible, wet down all dirty cages before dumping soiled bedding.
- Dump all soiled bedding in a suitable Bedding Disposal Workstation.
- Dispose of soiled bedding and other laboratory waste according to your facility's guidelines.

MAINTENANCE

- Workstations and safety cabinets should be re-certified by trained technicians at least once every 12 months.



WORLD CLASS. WORLDWIDE.