

# invent-uk

ANSI/ASHRAE 110-1995

Flow Containment & Instrumentation

## CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH ASHRAE 110-1995

CERTIFICATE & REPORT NO: INV/ASHRAE110/608

DATE: 30th November 2012



**Fume Hood Manufacturer:**

**Esco Micro Pte Ltd**  
**21 Changi South Street 1**  
**Singapore 486 777**

**Fume Hood Type:**

**1.8m wide EFD-6B1, EFD-6B2, EFD-6B3 bench-type**

**Sash Opening Width: 1600mm**

**Sash Opening Height: 440mm (100%)**

**Average Face Velocity: 0.51m/s (at 100%)**

**This is to certify that the fume hood described above has been type-tested in accordance with ASHRAE 110-1995 and resulted in performance characteristics given in the corresponding test report.**

**Tested and Certified by: Dr A F Bicen**

**Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR**

**Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com**

A handwritten signature in blue ink that reads 'A.F. Bicen'.