

**CERTIFICATE OF TYPE TESTING
IN ACCORDANCE WITH ASHRAE 110-1995**

**CERTIFICATE & REPORT NO: INV/ASHRAE110/365
DATE: 6th September 2006**



Fume Hood Manufacturer:

**Esco Micro Pte Ltd
21 Changi South Street 1
Singapore 486 777**

Fume Hood Type:

**1.2m wide EFB-4A1 bench-type
Sash Opening Width: 1050mm
Sash Opening Height: 460mm (100%)
Average Face Velocity: 0.50m/s (at 100%)**

This is to certify that the fume hood described above has been type-tested in accordance with ASHRAE 110-1995 and resulted in velocity and containment performance given in the corresponding test report.

Certified by: Dr A F Bicen

**Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR
Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com**

A handwritten signature in blue ink that reads 'A.F. Bicen'.