




Isolate the microscope from vibrations generated from the external environment



Note:
AVT-1 is designed to be easy to use and almost maintenance free

Technical Specifications	AVT-1
Overall dimensions (W x D x H)	1200 x 800 x 800 mm (47.25" x 31.5" x 31.5")
Weight	70 kg (154 lbs)
Material	Material Powder painted mild steel, Stainless steel and glass
Float dimensions (W x D x H)	540 x 340 mm (21.26" x 13.39")
Resonance frequency	5,5-50 Hz
Recommended Load weight	26~32 kg



 Designed in Denmark
 Made in the E.U.

900109_AVT_Equipment_AVT1_Files_V02_A4_Mar2014
Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials. Esco reserves the right to change specifications and prices without notice. All trademarks and logos are the property of Esco and the respective companies.