

Key Features

- ✓ Constantly validate up to 6 x CO₂ / O₂ incubators
- ✓ Controllable flow rate
- ✓ Monitor up to 6 x PT1000 sensors
- ✓ 6 ports for sequential gas samples
- ✓ Gas feedback returns sampled gas to incubator or exhaust

As validation unit, Miri GA can provide

- ✓ Parameter data logging
- ✓ Audible and visual alarms
- ✓ Adjustable alarm thresholds

Connected to a PC, the logging software will document parameters and optionally send e-mail alarms if set up.



 Designed in Denmark

 Made in the E.U.

Certified  EN 61010-1:2010

Specifications

Input ports	6 x PT1000 ports for temperature monitoring 6 x gas sampling ports
Output ports	1 x gas feedback port 1 x USB port
Shipping dimensions and weight	440mm x 430mm x 240mm (17.3" x 16.9" x 9.4"), 15kg (33.1lbs)

Ordering information ***

MRI-GA	Miri GA CO ₂ / O ₂ / Temp validation Unit, 115/230V, 50/60Hz
--------	--

Accessories

MRA-1101	1 pc PT1000 cable
MRA-1102	Set of 6pcs PT1000 cables
MRA-1103	1 pc Gas connection tube
MRA-1104	Set of 6pcs Gas connection tubes

*** Includes data logger software, 1 pc PT1000 cable, 1 pc Gas connection tube, 1 pc Gas feedback tube