

invent-uk

ANSI/ASHRAE 110-1995

Flow Containment & Instrumentation

CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH ASHRAE 110 1995

CERTIFICATE & REPORT NO: INV/ASHRAE110/325

DATE: 15th August 2005



Fume Hood Manufacturer:

Esco Micro Pte Ltd
21 Changi South Street 1
Singapore 486 777

Fume Hood Type:

1.34m wide ADC-4A1 bench-type
ductless recirculatory filtration
Sash Opening Width: 1270mm
Sash Opening Height: 200mm
Average Face Velocity: 0.53m/s

This is to certify that the fume hood described above has been type tested in accordance with ASHRAE 110 1995 and resulted in velocity and containment performance given in the corresponding test report.

Certified by: Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR
Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

Handwritten signature of A. F. Bicen in blue ink.