

invent-uk

BS 7258

Flow Containment & Instrumentation

CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH BS 7258

CERTIFICATE & REPORT NO: INV/BS7258/324

DATE: 15th August 2005



Fume Cupboard Manufacturer:

Esco Micro Pte Ltd
21 Changi South Street 1
Singapore 486 777

Fume Cupboard Type:

1.34m wide ADC-4A1 bench-type
recirculatory filtration

Sash Opening Width: 1270mm

Sash Opening Height: 200mm

Average Face Velocity: 0.54m/s

This is to certify that the fume cupboard described above has been type tested in accordance with BS 7258 Part 1 & 4 and resulted in performance characteristics given in the corresponding test report.

Certified by: Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR

Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A handwritten signature in blue ink that reads 'A.F. Bicen'.