

invent-uk

BS 7989

Flow Containment & Instrumentation

**CERTIFICATE OF TYPE TESTING
IN ACCORDANCE WITH BS 7989
CERTIFICATE & REPORT NO: INV/BS7989/327
DATE: 15th August 2005**



Filter & Fume Cupboard Manufacturer:

**Esco Micro Pte Ltd
21 Changi South Street 1
Singapore 486 777**

Filter Type:

2 x activated carbon - Indocarb-D612

Width: 570mm

Depth: 510mm

Height: 60mm

Fume Cupboard Type:

1.34m wide ADC-4A1 bench-type recirculatory filtration

Sash Opening Width: 1270mm

Sash Opening Height: 200mm

Average Face Velocity: 0.54m/s

This is to certify that the gaseous phase filter described above has been type tested in the fume cupboard described in accordance with BS 7989 and resulted in performance characteristics given in the corresponding test report.

Certified by: Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR

Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A handwritten signature in blue ink that reads 'A.F. Bicen'.